



TTAF[®]
medikal

FL Forehead
Infrared Thermometer



For a healthier world
we research, develop and produce!



TTAF[®]

www.ttafmedikal.com

Kavaklı Mh. İstanbul Cd. No:21
Beylikdüzü -İstanbul /TÜRKİYE

Tel : +90 212 855 82 20 (Pbx)

Fax : +90 212 855 82 21

©     /ttafelectronic

FL Forehead Infrared Thermometer

Align the forehead with a distance of about 1cm and press the mode button
Get the measurement results immediately


Measure ear canal temperature or forehead temperature dual mode

Forehead Gun Mode Non-contact temperature measurement

Won the "2019 Golden Dot Design Award Badge



Technical Specifications



Product Dimensions : 4.2cm*13.5cm*5.5cm
Net Weight : 47g
Battery : 2 Adet AA Batary
Measurement Range : 34.0 C (93.2F)
-42.0 C (107.6 F)